

Kinderen Stem Learning Insect Present Observation Kindermicroscop Educatief Microscope 1000X WiFi Bio Bio Microscop

Productomschrijving



* M1 is een intelligente microscoop met meerdere interfaces: WiFi en Type-C. Het kan ook worden verbonden met telefoon, pad, laptop en pc, met afbeeldingen direct zonder ogen te gebruiken. Het is gemakkelijk in foto en video via de app en sla de afbeeldingen op. * M1 heeft geen scherpe en duidelijke afbeeldingen, maar ook handig om te nemen vanwege de dedalaniteit. Het heeft de gereflecteerde LED-lichtbron en de basis heeft 7-kleuren verstelbare lichtbronnen zonder kracht, van violet tot wit en u kunt vrijelijk de uitgezonden lichte kleur inschakelen. Elimineert de fase van kunstmatige verven voor transparante en doorschijnende secties. * Upload via de Type-C-interface, laden zonder het netsnoer. Ondersteun het gebruik van de lading van de lading PAL direct, de 2700 mAh hoge dichtheid lithiumbatterij, de basis kan worden gebruikt om tot 10 uur mee te gaan en de monomeermicroscoop kan afzonderlijk tot 2 uur worden gebruikt. * Deze intelligente microscoop kan op verschillende gebieden worden gebruikt: PCB-inspectie, horloge en telefoonreparatie, juwelenidentificatie, leren en demonstreren, textieltests, huiddetectie, afdrukinspectie, geldbesturing enz.

Specificaties

Garantie	1 jaar
Aangepaste ondersteuning	OEM, ODM.
Theorie	Videomicroscoop

Drawtube.	Monoculair
Totale vergroting:	1x ~ 1000x.
Materiaal:	ABS, Aluminiumlegering, PMMA
Microscope-batterij:	700mAh.
Basic Battery Stand:	2000mAh.
Koppel:	WiFi, Type-C, Venster, Android
Reflected Lighting:	LED-lichtbron.
Basic Lighting:	7 kleuren lichtbron verzonden
Lichtbron:	8 LED-verlichting
Kenmerkend:	Draadloze microscoop.
Certificering:	REACH, ROHS, CE, FCC, KC, PSE
Bruto gewicht	2.7kgs.
Dimensie	25 * 25.5 * 21 cm





Explore the Micro World Around You



Bee wings



Ant



Bee leg



Longitudinal section
of corn seed



Slice ovary
of lily



Longitudinal section
of pumpkin stalks









How to Operate on



For Win10

Unpack and pick out microscope, connect the microscope and host computer via USB cable, search Camera in Windows for find it and click on.

(Notice: please see user manual for details)

For Win7 / Win8 / Mac

Step 1. Install software "HiView"

Step 2. Connect the microscope to PC via the USB cable, double-click to open the software

Step 3. Click "Device" in the menu bar, and then click "Open", select the microscope name "*****", and click "Open" to enter the preview screen

How to Use on Cellphone



Download the APP.



Long press the power button to turn on the microscope.



Turn on the phone's WIFI setting, find the name: "wifi_camera_wifi30_xxxx" and connect it.




Then you can click the APP and use



How to Use on Cellphone

 (Notice: please see user manual for details)



Step 1. Download the APP.



Step 2. Long press the power button to turn on the microscope.



Step 3. Turn on the phone's WIFI setting, find the name: "wifi_camera_wifi30_xxx" and connect it.



Step 4. Then you can click the APP

Wide Compatibility

Wi-Fi and USB 2 in 1 connection



iPhone



Android



Tablet



Mac




Windows



Wi-Fi



USB cable





Verpakking en levering

