Ear Cleaning Endoscope



Specifications

- 1. 3 in 1 USB electronic Otoscope Borescope Ear Camera: The product is a high-definition visual ear cleaning tool ,convenient to use on OTG Android device, tablet,Windows PCs ,Windows XP/7/8/10(Not support iPhone)
- 2. Interesting Household Tool: It can tell the ear canal health status and not only can check the ear canal, but also can check other body parts, such as scalp, pores, comedones, ear canal, oral cavity, nasal cavity and so on.
- 3. Image Capture and Video Recording:It can take pictures and videos on the software interface,and you can browse pictures and videos.The best focal length is about: 1.5cm
- 4. Waterproof, dust-proof, 6 pcs high brightness LED light.
- 5. One-piece design, easy to carry, the visible inside of the ear spoon change your mind more than a little. Just connect the USB with your phone or PC, then put it onto the head of camera, download the APP/Software then see the live image.

Product Pictures:



Ear Endoscope Installation Methods



1. Find the locating hole at the front end of the camera pen



2.Insert the bent part at the end of the metal ear spoon into the locating hole



3. Push the fixing sheath to the bent part of the ear scoop. Don't push it to the end, just stop when it's fixed.





Approval Certification:







Why Choose US:

Our goal is to provide customers products of high quality, competitive price, and superior services.

- 1.Established in 2005, about 10 years of history;
- 2.Own main factories, supply diversification products and One stop service
- 3. Certificates of ISO9001, UL, CE, EMC, CQC, RoHS, Reach, etc.
- 4. Six years Gold member in Alibaba, Aliexpress etc
- 5. Production capability more than 100 thousands pcs/ month;
- 6.Monthly export more than 50 times delivery;
- 7.Exported to more than 20 countries and regions;
- 8. We provide an excellent combination between prices, quality & respond, delivery, after-sales!

